



## ES256H GENERAL DISINFECTANT CONCENTRATE

- Phosphate-free formulation
- Neutral pH will not dull finishes with repeated use

	ES256
APPEARANCE	Clear blue liquid
ODOR	Lemon
SOLUBILITY	Soluble
RELATIVE DENSITY	1.02
pH-VALUE	7-8
SHELF LIFE	3 Years

SIZES	
ES256-CS	4.73 L (1.25 U.S. gallon)
ES256-4CS2	1.89 L (0.50 U.S. gallon)

### ACTIVE INGREDIENTS GUARANTEE

Octyl decyl dimethyl ammonium chloride.....	6.510%
Diocetyl dimethyl ammonium chloride.....	3.255%
Didecyl dimethyl ammonium chloride.....	3.255%
Alkyl (50% C14, 40% C12, 10% C16) dimethyl benzyl ammonium chloride.....	8.680%

### DILUTIONS

USE	RATIO	PER GALLON OF WATER	PER LITRE OF WATER
Disinfectant/ Sanitizer/Cleaner	1:256		

**CAUTION:** Refer to label for complete instructions, dilutions, safety requirements and kill claims

### BENEFITS

- One step neutral disinfectant that is effective against a broad spectrum of bacteria, is virucidal\* (including HIV-1, HIV-2, HBV and HCV) and inhibits the growth of mold and mildew and their odors when used as directed
- Cleans, shines, deodorizes and disinfects all hard non-porous surfaces listed on the label
- Neutral pH multi surface cleaner, deodorizer and disinfectant

### USE

Use on floors, walls, toilet bowls, urinals and other hard nonporous surfaces such as sinks, tables and counter tops made of Stainless Steel, Metal, Glazed Porcelain. Other areas of use include:

- Farms
- Florists
- Hospitals/Nursing Homes
- Hotels
- Kitchens
- Medical, dental and veterinary offices
- Schools

Can be used with mop and bucket, trigger sprayer, sponge soaking or by auto scrubber.

### DISPOSAL

Dispose of used/unwanted product according to state and local regulations. Offer container for recycling.

### PERSONAL PROTECTION

Please refer to product label and SDS